

Datasheet

MORF4L2 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009643-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant MORF4L2.

Immunogen: MORF4L2 (NP_036418.1, 1 a.a. ~ 288 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MSSRKQGSQPRGQQSAEEENFKKPTRSNMQRSKMR
GASSGKKTAGPQQKNLEPALPGRWGGRSAENPPSG
SVRKRTRKNQKTPGNGDGGSTSEAPQPPRKKRARAD
PTVESEEAFFKNRMEVKVKIPEELKPWLVEDWDLVTRQ
KQLFQLPAKKNVDAILEEYANCKKSQGNVDNKEYAVN
EVLVAGIKEYFNVMLGTQLLYKFERPQYAEILLAPDAP
MSQVYGAPHLRLFRIGAMLAYTPLDEKSLALLLGYL
HDFLKYLAKNASLFTASDYKVASAEYHRKAL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9643

Gene Symbol: MORF4L2

Gene Alias: KIAA0026, MORFL2, MRGX