

Datasheet

MTHFS monoclonal antibody (M01), clone 2C12

Catalog Number: H00010588-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant MTHFS.

Clone Name: 2C12

Immunogen: MTHFS (NP_006432, 104 a.a. ~ 203 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

LPKTSWNIPQPGEGDVREEALSTGGDLIFMPGLGFD
KHGNRLGRGKGYDAYLKRCLQHGEVKPYTLALAFKE
QICLQVPVNENDMKVDEVLYEDSSTA

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10588

Gene Symbol: MTHFS

Gene Alias: FLJ30410, HsT19268