



Certificate of Conformity

VWR CLEANROOM WIPES, NON WOVEN, SPEC-WIPE® 3

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Description	Sterility Level	L×W (mm)	Packed	Pk
115-0597	Spec-Wipe® 3	Non-sterile	229×229	300/bag	300
115-0598	Spec-Wipe® 3	Non-sterile	305×305	150/bag	150

Material: Durable, non woven wipes, made from 46% polyester and 54% cellulose.

Attribute; (units)		Value **	Test Method
Basis weight; (g/m ²)		52	
Absorbency in water	Intrinsic; (mL/g)	5.43	IENT-RP-CC004.2, Sec. 7.1
	Extrinsic; (mL/m ²)	283	IENT-RP-CC004.2, Sec. 7.1
Sorptive rate; (seconds)		1	
Extractables	In deionized water; (g/m ²)	0.025	IENT-RP-CC004.2, Sec. 6.1.2
	In isopropyl alcohol; (g/m ²)	0.020	
Specific ions	Sodium; (ppm)	8.3	IENT-RP-CC004.2, Sec. 6.2.2
	Chloride (ppm)	9.4	
Particle generation Po ≥ 0.5µm; (x10 ⁶ /m ²)		58.2	IENT-RP-CC004.2, Sec. 5.1
Fibers > 100µm; (x10 ³ /m ²)		44	IENT-RP-CC004.2, Sec. 5.2

** ND = None detected; levels are below detection limit of test equipment



Certificate of Conformity

Other Characteristics:

Hydroentangled: Yields a durable wiper which produces less particulation than some other nonwoven methods
Hydroentangled without the use of chemicals: Reduces ionic contamination
Lightweight polyester/cellulose combination: Economical and sorbent
Creped fabric: Provides improved softness and bulk

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **January 2018**