



Pharma. Diluent USP

AX021203

For the dilution of samples in the analysis of non-sterile products for microbial contaminants

Presentation	Packaging Details	Shelf Life (Months)	Storage
Flasks 300 ml	12 x 300 ml	12 month	15-25°C

COMPOSITION (g/l)	USING
sodium chloride8.5	Depends on the purpose for which the media is used.
disodium hydrogen phosphate	
dodecahydrate0.29	
sodium dihydrogen phosphate0.16	

QUALITY CONTROL

Physical/Chemical Control

pH: 7.2 ± 0.2 at 25 °C

Microbiological Control

Test strains	Enumeration at T0 and T45 at 20-25°C
Escherichia coli ATCC 25922	±30% / T0
Staphylococcus aureus ATCC 6538P	±30% / T0

According to
USP

BIBLIOGRAPHY